


Profilyzer

Determination of sugars and alcohols

SunData wine system

Method : Zucker
 Method table : zuck0709.tbl
 Sample : Std blau
 Sample number : 53
 Concentration factor : 1.0000
 Standard : c:\sundata\data\Z0709036.RPT

SunChrom  [®]
 ...mit uns auf der Sonnenseite der Chromatographie



Peak	Retention time	Substance	Peak Height	Amount [g/l]	Vol%
1	0.40		0.589		
2	1.55	Glucose	60.553	32.292	
3	2.65	Fructose	59.719	33.864	
4	4.28	Glycerol	17.344	9.085	
5	5.42	Ethanol	69.857	81.244	10.294

Total sugar : 66.157 g/l
 Total alcohol : 112.338 g/l (14.233 Vol%)
 Grape juice addition: 21.888 Vol%